


Masterfilter



Filters for the process industry

Masterfilter



Masterfilter was created in 2010. Competence in processes and utilization are the strengths of Masterfilter.

We bring solutions for all kind of critical applications.

Masterfilter is a filtration company that focuses on high end solutions in PES & PTFE membrane filters and various pleated filter cartridges. Masterfilter provides professional technical and scientific assistance to filtration, separation and purification solutions. Our goal is to deliver innovative quality products on time with a high cost saving potential.

The philosophy of Masterfilter

- Experience
- Market and customer oriented
- Flexible
- Team spirit
- High qualified engineers
- Stable and independent
- Very competitive
- Find best fitting solution



Filters for the process industry

MF
Masterfilter



Wine & Cider

Activated Carbon Filter

ACF



- Ideal for chlorine reduction
- Excellent absorbency
- No release of active carbon
- Available in 2 grades
- FDA approved

Applications

- De-chlorination of process water and drinking water
- De-chlorination of drinking water
- De-chlorination of rinsing water
- De-colourisation

Melt blown depth filter

MBK



- Three layer construction
- Thermal bonding without binders and adhesives
- High dirt holding capacity and long service life
- Removal rating (μm) : 1,0 – 150 μm

Applications

- Pre-filtration for wine
- Pre-filtration for RO plants
- Beverage filtration

Multi-layer pleated polypropylene cartridge

APK



- Multi layers melt-blown polypropylene
- High dirt capacity
- 100 % polypropylene
- Absolute micron rating (μm) : 0,3 – 10 μm
- Retention of particles
- Reduction of bioburden from liquids

Applications

- Pre-filtration for wine, water and gases
- Clarification of wine

Nano-fiber polypropylene pleated filter cartridge

APF



- High flow rate
- High dirt capacity
- Back washable, extend service life
- Micron rating (μm) : 0,3 – 20 μm
- Absolute micron rating (μm) : 1,5 – 10 μm

Applications

- Trap filtration for wine
- Wine clarification
- Yeast removal from wine
- Trap filtration for kieselguhr free microfiltration

High density polypropylene pleated cartridge

PPK



- Graded pore-size distribution
- Back washable, extend service life
- 100 % polypropylene
- Micron rating (μm) : 0,3 – 20 μm

Applications

- Pre-filtration for wine, water and gases
- Trap filtration for wine

Hydrophilic polyethersulfone membrane

PES



- Easy wetting
- Low adsorption
- Repeatedly integrity testable
- Micron rating (μm) : 0,1 – 1,2 μm

Applications

- Bottled water
- Wine filtration
- Sterile wine filtration
- Wine



Sterile PTFE membrane gas filter

SLF-T



- Permanently hydrophobic PTFE membrane
- Suitable for all integrity test procedures
- Pleated surface for high flow rate and long life
- Removal rating (μm) : 0,2 μm

Applications

- Sterile filtration of air and CO_2
- Sterile tank venting
- Process gas filtration

Glass micro fibre cartridges

SLF-G



- Steam sterilisable
- High dirt holding capacity
- Pleated surface for high flow and low pressure drop
- Removal rating (μm) : 0,01, 0,3, 0,5 μm

Applications

- Sterile filtration of air and CO_2
- Sterile tank venting
- Process gas filtration
- Membrane pre-filter

Stainless steel sintered filter cartridges

SDF



- Available in 316 stainless steel
- Cleanable and re-usable
- Welded construction
- High temperature and corrosive resistant
- Removal rating (μm) : 1 – 100 μm

Applications

- Steam filtration (culinary)
- Aggressive chemicals and gases
- Viscous fluids
- Corrosive fluids

Housing range for liquids and gases

FG-F/G



- Manufactured in accordance with European regulations
- Supplied with a wide range of connections
- Mechanically / electro polishing
- Surface finish $R_a < 0,8 \mu\text{m}$
- Sanitary / industrial

Applications

- Liquid filtration
- Gas filtration
- Vent filtration

Filter integrity tester

IT-01B



- Integrity testing for all membrane filters
- Single and multiple test
- Liquid and gas filter
- Compact design

Applications

- Bubble point test
- Pressure decay test
- WIT (optional)

Masterfilter

Belgium

Masterfilter Belgium
Groddeveg 14 /A
B-3010 KESSEL-LO
Belgium

info@masterfilter.be
www.masterfilter.eu
Tel: +32 13 67 09 04
Fax: +32 13 67 09 01

Germany, Austria & Switzerland

Masterfilter GmbH
Kapellenstrasse 9
85622 Feldkirchen
Germany

info@masterfilter.com
www.masterfilter.com
Tel: +49 89 2885 8787

The Netherlands

Masterfilter Nederland
W.J. van der Hoofstraat 18
4194 WD Meteren
The Netherlands

info@masterfilter.nl
www.masterfilter.nl
Tel: +31 345 50 83 96

Your local distributor



Filters for the process industry

MF
Masterfilter